



# precision antibody™

- Jun Hayashi, Ph.D.
- President
- 9130 Red Branch Rd
- Columbia, MD 21045
- U.S.A.
- [jhayashi@precisionantibody.com](mailto:jhayashi@precisionantibody.com)
- +1-410-884-4100

# 60日で完全ヒトモノクロ抗体の作製



Trans **C**hromosomes



## 戰略的提携

**Precision Antibody**, a service division of A&G Pharmaceutical Inc is pleased to announce business alliance with **Trans Chromosomics**



### A&G Pharmaceutical Inc/Precision Antibody

- Theranostic antibody development against breast & lung cancer
- Global leader in antibody development



CEO, Ginette Serrero Ph.D.  
Established: 1, June 2000  
[www.agpharma.com](http://www.agpharma.com)



President, Jun Hayashi Ph.D.  
Established: 1, June 2000  
[www.precisionantibody.com](http://www.precisionantibody.com)



### Trans Chromosomics (TC)

- *Global Leader Of Artificial Chromosome Engineering Technology*
- *Innovative Platforms For Biopharmaceuticals, Cell/Gene Therapy, And Xenotransplantation*




President CEO, Mitsuo Oshimura Ph.D.  
Established : 17 December 2014  
[www.trans-chromo.com](http://www.trans-chromo.com)





## 概要

- 完全ヒト抗体産生マウス（TC-mAb™マウス）は、免疫抑制がなく、 Precision Antibody社（PA）の技術で短期間に、力価の高い抗原特異的抗体を作製可能
- PAは、TC-mAb™マウスから、約60日で完全ヒトモノクローナル抗体（mAb）を作製可能
- TC-mAb™マウスは100%ヒト抗体を発現
- 作製された完全ヒトmAbは、抗体遺伝子のヒト化といった改変が不要
- 作製された完全ヒトmAbは、正常なヒトIgGレパトアを保持
- 抗原特異的完全ヒトmAb産生ハイブリドーマ細胞を高効率で作製可能
- 多数のクローン・多様なエピトープを保持する抗体を作製可能



## 事例紹介

### 事例1

免疫: ~90KDa protein with ~75% homology with mouse ortholog

期間: 60 days from immunization to expansion of antigen-specific clones.

抗血清の力価（融合前）：

Serum Dilutions	1:1K	1:3K	1:10K	1:30K	1:100K	control
A <sub>620</sub>	3.40	2.90	1.66	0.64	0.10	0.04

ELISAによる抗原特異的抗体産生細胞数の評価:

antigen: 30 ng/well      1°Ab: 100 ng/well

#### OD # of Ag-positive clones (out of 1727 clones)

0.5-1.0	335
1.0-2.0	423
2.0->3.0	387
>3.0	183
<b>Total Positive</b>	<b>1,228</b>

~77% of clones were positive against the immunogen.

### 事例2

免疫: ~50KDa protein no homology to mouse proteins

期間: 60 days from immunization to expansion of antigen-specific clones.

抗血清の力価（融合前）：

Serum Dilutions	1:1K	1:3K	1:10K	1:30K	1:100K	control
A <sub>620</sub>	3.02	2.94	2.60	1.80	0.95	0.04

ELISAによる抗原特異的抗体産生細胞数の評価:


antigen: 50 ng/well      1°Ab: 100 ng/well

#### OD # of Ag-positive clones (out of 1920 clones)

0.5-1.0	358
1.0-2.0	205
2.0-3.0	276
>3.0	118
<b>Total Positive</b>	<b>957</b>

~50% of clones were positive against the immunogen.






## TC-mAb™マウスの治療用抗体開発に対する優位性


	Wild-type	Fully human IgG
<b>Animals</b>	Balb/c, ICR, etc.	TC-mAb mice
<b>Immune response</b>	Mouse IgG	Human IgG
<b>Antibody titer</b>	Excellent	Excellent
<b>IgG subclass</b>	mlgG1, 2a, 2b, and 2c	hlgG1, 2, 3, 4
<b>Humanization</b>	The entire region	Not necessary
<b>Ag-specific mAb generation</b>	Normal	High ratio
<b>Size of spleen</b>	Enlarged	Normal
<b>B cell development</b>	Normal	Expanded
<b>Antigen-specific B cells</b>	Normal	Expanded
<b>Ig-gene Repertoire</b>	Normal mouse	Normal human
<b>Somatic hypermutation</b>	Yes	Yes
<b>Length of CDR3H</b>	Mouse (short)	Human (long)

### Summary


- 正常なヒトIgGレパートアを保持
- 完全なヒト抗体かつ、親和性が成熟しているため、その後の改変が不要



## TC-mAb™マウスにおけるヒト抗体遺伝子のレパートア解析 I





Heat map representation of VH and V<sub>κ</sub> -gene usage




## TC-mAb™マウスにおけるヒト抗体遺伝子のレパートア解析Ⅱ

## VH gene usage (IgM&amp;IgG)


VH gene usage (IgM &amp; IgG)

V $\kappa$  gene usageV $\kappa$  gene usage




## TC-mAb™マウスにおけるヒト抗体遺伝子のレパートア解析Ⅲ

CDR3H Length



CDR3H Amino acids usage



Both amino acid length and composition in CDR3 region of the antibody heavy chain variable region (VH) were almost the same with human PBMC.

It was speculated that antibodies similar to human individuals can be produced from TC-mAb mice.